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70*90mm

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JTF Just To Future

J40 ELECTRIC TOOTHBRUSH

USER MANUAL

Product Specification	
Product name	Electric Toothbrush
Model No.	J40
Output	1.8W
Rated input	3.7V
Battery last time	45 days of use (when fully charged)
Waterproof Level	IPX7
Package include	
Brush Handle x1, Brush Head x6, Tongue Scraper x1, Travel Bag x1, USB DC Charging Cable x1, User Manual x1	

Product Description

- ①Brush Head
- ②Brush Handle
- ③Power On/Off button
- ④Brushing mode
- ⑤Battery Indicator
- ⑥Charging port

Features Brushing modes

J40 electric toothbrush is equipped with 5 different modes (shown as picture)
Mode selection:
--Click the power button to turn on the toothbrush.
--Click the button again to cycle between modes. Choose the mode you need.
--If there are no more operations, the toothbrush will automatically shut down after 2 minutes. You can also long press the button for 5 seconds to turn off.
--The toothbrush will automatically remember and stay in the mode the next time you turn it on.

Smartimer

The toothbrush will automatically shut down after 2 minutes, indicates that your brushing cycle is complete. Dental professionals recommend brushing no less than 2 minutes twice a day.

Quadpacer

The Quadpacer is an interval that has a short pause every 30 seconds to remind you to brush the 4 sections of your mouth evenly and thoroughly.

Charge

- ①Turn on the toothbrush before using it for the first time. If you feel it is not powerful enough, recharge it.
- ②Charging for 2 hours the toothbrush can be used for about 45 days.
- ③If the toothbrush has not been charged for a long time, please wait 6 to 10 minutes for the indicator light to turn on.
- ④If the toothbrush has not been used for more than 6 months, charge it for at least 12 hours.

Brush head

- ① Rinse the brush head and bristles after each use.
- ② Remove the brush head from the handle and rinse the brush head

connection with warm water at least once a week. Do not clean with hot water above 40°C or acidic cleaning solutions.
③Replace the brush heads every 3 months for optimal results. Use only J40 brush heads compatible with JTF J40 brush handle.

Battery Indicator

- When the battery is low: the battery indicator will flash red to remind you to charge.
- When the battery is less than 3%: the battery indicator will flash red for a few times and then automatically shut down.
- During charging, the red light flashes. When charging is completed, the red light stays on.

Fault Diagnosis

Problem	Possible Causes	Counter-measure
Reduction of vibration intensity	① The brush head puts too much pressure on the teeth ② The brush head is too close to the handle	① Place the brush head lightly against your teeth. ② Remove the brush head from the handle and re-set it on the output shaft, keeping a distance between the toothbrush head and the output shaft handle
Even after charging, the toothbrush can only be used for a short time.	① Insufficient charging time ② This product has not been used for more than 6 months ③ The battery has reached the end of its service life.	① Charge until the indicator light stays red ② Charge the toothbrush continuously for 12 hours ③ Contact after-sales service
Toothbrush doesn't work	① Insufficient battery ② The battery has reached the end of its service life	① Charge the toothbrush for 2 hours ② Contact after-sales service

If you have any questions, please email us at support@jtforal.com

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