

TOARNA conținutul unui stick în 200 ml apă rece	AMESTECA până la dizolvare completă	CONSUMA imediat după prepararea soluției
---	---	--

se agită până la omogenizare și se consumă imediat după prepararea soluției.

MOD DE UTILIZARE: Adultii 1-2 sticle unită. Conținutul unei sticle se diluează în aproximativ 200 ml apă rece.

HYDRONIO UBE - Stick-uri
SUPPLEMENT ALIMENTAR CU INDUCITOR
28 stick-uri (cu aroma de aine)

INFORMATION NUTRITIONALE:

Ingrediente	Cantharol/500k	VRB+ %	2 stocuri
Solui	160 mg	0	-
Polisid (material solidat în cantitate necunoscută de mg) (15/18/19)	150 mg	0	16%
Magneziu	175 mg	33.3 %	66.7 %
Polisid (Clorid ascorbic)	50 mg	42.5 %	17.5 %
Polisid (Clorid sub formă de acid acetic)	5 mg	3%	67%
Zinc	3 mg	30%	60%
Polisid (Clorid paracetamol)	2 mg	33.3 %	66.7 %
Polisid (Clorid paracetamol)	0.3 mg	21.4 %	42.8 %
Polisid (Clorid paracetamol)	5 mg/200 U.L.	100%	200%
Polisid (Clorid paracetamol)	5 mg	200%	400%

Yates, 2004; Mays et al., 2004; Pech et al., 2004; Regalado et al., 2004; O'Connell et al., 2005; Johnson et al., 2006).

HYDRONI**Q**
UBE

ELECTROLITŢI

+ DEXTROZĂ

SUPLIMENT ALIMENTAR CU ÎNDULCITOR - 28 STICKS x 7.5 g

HYDRONID		ELECTROLITI		MINERALE		VITAMINE	
----------	---	-------------	---	----------	---	----------	---

HYDRONIQ

ELECTROLITIC

MINERALE

VITAMINE

 Laboratoarele
Remedia

The 5 women interviewed were 23–30 years old and had been

inhibition against the growth of *S. aureus*.

Produk resmi dari: CASIO COMPUTER CO., LTD. 2024

Green tube weight: 275 g (20 stick-out (cattle))

A se constata, de probetura, tratente de integritate laudat:

5